

Acu-Form® Penile Prosthesis



Malleable Penile Prosthesis

- Superior positionability, virtually eliminates springback
- Excellent distal shaft column strength helps prevent buckling
- No springs, cables or moving internal parts to compromise reliability
- Easy assembly, no special tools required
- Trimmable: custom fit to each individual's total corporal length
- A minimal inventory of three prostheses can cover all your implant needs



Lifetime Replacement Policy


Acu-Form® Penile Prosthesis

Catalog Number	Diameter Size	Description
9509S	9.5mm	One pair penile prostheses and three pairs tail caps with connectors (standard, +.5 cm and +1 cm). The implanted length can vary from 14 cm to 23 cm (when the +1 cm tail cap is used).
9511S	11mm	One pair penile prostheses and three pairs tail caps with connectors (standard, +.5 cm and +1 cm). The implanted length can vary from 16 cm to 25 cm (when the +1 cm tail cap is used).
9513S	13mm	One pair penile prostheses and three pairs tail caps with connectors (standard, +.5 cm and +1 cm). The implanted length can vary from 18 cm to 27 cm (when the +1 cm tail cap is used).
9320S	9.5mm	Malleable Sizer. Consists of one 9.5 mm diameter cylinder marked in .5 cm increments from 2 cm - 20 cm. Available nonsterile.

Prosthesis

Trim marks indicate final length with standard tail cap.

Caution: Do not trim the prosthesis past the most distal mark on the shaft.



Sizer

Note: The reverse side is calibrated in the opposite direction.



Tail Cap



Measuring: The centimeter marks indicate the final length of the prosthesis after the standard tail cap is firmly attached to the proximal end (e.g., trimming at the 18 cm mark, plus standard tail cap, equals 18 cm total length).

Returns

Authorization must be obtained from Mentor prior to the return of device. Merchandise may be returned for exchange or credit within 30 days of invoicing. Returns may be subject to a restocking charge.

Ordering Information

To order, please contact the Mentor Customer Service Department at Mentor, 201 Mentor Drive, Santa Barbara, California 93111, U.S.A. Telephone (805) 879-6000

Note: The Acu-Form® Penile Prosthesis is supplied sterile. The device should not be resterilized.



LIFETIME PRODUCT REPLACEMENT POLICY FOR MENTOR[®] PENILE IMPLANTS

Alpha I[®] Inflatable Penile Prosthesis
Mark II[®] Inflatable Penile Prosthesis
Acu-Form[®] Penile Prosthesis
Malleable Penile Prosthesis

Mentor will replace the Mentor Penile Prostheses noted above or any component thereof for any reason during the lifetime of the patient-recipient of a Mentor Penile Prosthesis upon request. Requests for replacement product or component thereof should be made by the patient's surgeon to Mentor at the toll-free number listed below, at which time a Return Authorization Number will be provided. Mentor will then provide a replacement Mentor Penile Prosthesis or component thereof without charge. Any Mentor Penile Prosthesis or component so replaced should be returned to Mentor and shall then become the property of Mentor.

This replacement policy is limited to replacement of the Mentor Penile Prostheses or components and Mentor shall not be responsible for any incidental or consequential damages of any kind, including without limitation, medical, hospital, surgical and hospital costs or the reimbursement for the cost of penile implants not manufactured by Mentor. This product replacement policy shall be effective provided Mentor manufactures penile prostheses at the time the request for replacement is made. This policy is not a warranty or guaranty. Mentor excludes all warranties, whether express or implied by operation of law or otherwise, including without limitation, any implied warranties of merchantability or fitness. Mentor shall have no responsibility for any statements or representations, whether oral or written made by any other person or entity.

Mentor makes no representation as to the life expectancy of the penile prosthesis. There are conditions which may cause these mechanical devices to wear out or malfunction and require replacement.



5425 Hollister Avenue • Santa Barbara CA 93111
Telephone (805) 681-6000 • Toll Free (800) 235-5731

© May 1998 Mentor
0294/US05
9805011



Corporate Offices
201 Mentor Drive
Santa Barbara, California 93111
Tel: (805) 879-6000
Fax: (805) 967-7108
www.mentorcorp.com

**To order product please call:
Customer Service (800) 235-5731**

International Sales Offices

Mentor Medical Systems Australia, Pty. Ltd.
Unit 2/935 Station Street
Box Hill North
Victoria 3129
Australia
TEL: (61) 3 9890-9444
FAX: (61) 3 9890-9777

Mentor Deutschland GmbH
Ludwigsforum, Ludwigstrasse 45/Haus C
D-85399 Hallbergmoos
Germany
TEL: (49) 8116-00500
FAX: (49) 8116-005011

Mentor Medical Systems Canada, Inc.
1333 Boundary Road, 10
Oshawa, Ontario, L1J 6Z7
Canada
TEL: (905) 725-7763
Toll Free: (800) 668-6069
FAX: (905) 725-7340

Mentor Benelux B.V.
Pampuslaan 1
1382 JM Weesp
The Netherlands
TEL: (31) 294-491550
FAX: (31) 294-491555

Mentor Medical Systems Ltd., U.K.
The Woolpack, Church Street
Wantage, Oxon, OX12 8BL
United Kingdom
TEL: (44) 1235-768758
FAX: (44) 1235-768780

Mentor Medical Systems France, S.A.
Centre D'Affaire La Boursidiere
Le Plessis Robinson
92357 Cedex
France
TEL: (33) 1 47 21 33 66
FAX: (33) 1 47 29 99 83

Mentor Medical Systems Iberica, S.L.
C/Rodriguez Marin
88-Bajo DRCHA
28016 Madrid
Spain
TEL: (34) 91-562-2700
FAX: (34) 91-563-3951

Mentor Medical Italia, S.r.l.
Via Olmetto 7
20123 Milano, Italy
TEL: (39) 02-880-7761
FAX: (39) 02-7209-3710